

TABLE OF CONTENTS

1	Overview	1
1.1	date	1
1.2	motionTracker	1
1.3	parameterFile	1
1.4	rigCodeVersion	1
1.5	singleTrialTag	1
1.6	subjectName	1
1.7	tag	1
1.8	translationDate	1
1.9	trialName	1
1.10	trialNumber	1
1.11	trialStatus	1
2	Definitions	2
2.1	algTargTypeDef	2
2.2	AllPossibleAnalogTargets	2
2.2.1	dataSource	2
2.2.2	dataSourceType	2
2.2.3	linkedTargets	2
2.2.4	names	2
2.2.5	tolerances	2
2.2.6	windows	2
2.3	allPossibleIntertrialIntervals	2
2.4	AllPossibleIntervals	2
2.4.1	modes	3
2.4.2	names	3
2.4.3	primary	3
2.4.4	values	3
2.5	AllPossibleMarkers	3
2.5.1	channelNumber	3
2.5.2	dataSource	3
2.5.3	dataSourceType	3
2.5.4	names	3
2.6	AllPossibleStateTables	3
2.6.1	globalTimeout	4
2.6.2	globalTimeoutFailureStep	4
2.6.3	names	4
2.6.4	probabilities	4
2.7	AllPossibleTargets	4
2.7.1	formulas	4
2.7.2	modes	4
2.7.3	names	4
2.7.4	primary	4
2.7.5	selectionModes	4
2.7.6	shape	4
2.7.7	window	5
2.7.8	xValues	5
2.7.9	yValues	5
2.7.10	zValues	5

2.8	AllPossibleTubeFilepaths	5
2.9	AllPossibleTubeGroups	5
2.10	AllPossibleTubeLengths	5
2.11	AllPossibleTubeNames	5
2.12	AllPossibleTubePathLengths	5
2.13	AllPossibleTubeRenderTypes	5
2.14	allPossibleVibrotactorDefinitions	5
2.15	allPossibleVibrotactorDimensions	5
2.16	allPossibleWaveformShapes	6
2.17	analogChannelNames	6
2.18	autoMonkeyProfileDef	6
2.19	gazeCheckDef	6
2.20	handCheckDef	6
2.21	intervalModeDef	6
2.22	MarkerSpeedProfiles	6
	2.22.1 finalValues	6
	2.22.2 initialValues	6
	2.22.3 names	6
	2.22.4 relativeSpeed	6
	2.22.5 types	6
2.23	motorModeDef	7
2.24	rawDecodeColumnNames	7
2.25	snippetInfoRowNames	7
2.26	targetModeDef	7
2.27	targetSelectionModeDef	7
2.28	targetShapeDef	7
2.29	targetTypeDef	7
2.30	targetWindowShapeDef	7
2.31	ttlCheckDef	7
2.32	tubeCheckDef	7
2.33	vrPhasespaceTransformation	7
	2.33.1 rotation	7
	2.33.2 scaling	7
	2.33.3 translation	7
3	Parameters	8
3.1	eyeTrackerCalibrationMatrix	8
3.2	ForceParameters	8
	3.2.1 initialLambda	8
	3.2.2 inputMethod	8
	3.2.3 kicms	8
	3.2.4 kv	8
	3.2.5 kvib	8
	3.2.6 refValue	8
	3.2.7 slowfastRate	8
	3.2.8 systemName	8
	3.2.9 systemOrder	9
	3.2.10 visCursorOffsets	9
3.3	LockZParameters	9
	3.3.1 lockZ	9
	3.3.2 lockZCoordinate	9

3.4	markerJitter	9
3.4.1	A	9
3.4.2	controlJitter	9
3.4.3	S	9
3.4.4	visibleJitter	9
3.5	MotorParameters	9
3.5.1	amplitudeRanges	9
3.5.2	carrierFrequencyRanges	9
3.5.3	distanceRanges	10
3.5.4	dutyCycleRanges	10
3.5.5	forceToZeroIfOutsideRange	10
3.5.6	ignoreIfOutsideRange	10
3.5.7	modulationFrequencyRanges	10
3.5.8	vibrotactorDimensions	10
3.5.9	vibrotactorsNumber	10
3.5.10	vibrotactorTypes	10
3.5.11	waveformShapes	10
3.6	Pseudorandomness	10
3.6.1	pseudorandom	10
3.6.2	pseudorandomReset	10
3.7	StateTable	10
3.7.1	Analog	11
3.7.2	applyGlobalTimeout	11
3.7.3	autoMonkeySpeedProfileName	11
3.7.4	Display	11
3.7.5	Gaze	11
3.7.6	Hand	12
3.7.7	Interval	12
3.7.8	Motor	12
3.7.9	passTable	13
3.7.10	resetGlobalTimeout	13
3.7.11	stateName	13
3.7.12	StateTargets	13
3.7.13	togglePhotodiodeTarget	13
3.7.14	TTL	13
3.7.15	tubeColor	14
3.7.16	tubeMarker	14
3.7.17	tubeSection	14
3.7.18	tubeTexture	14
3.7.19	tubeWindow	14
3.7.20	tubeWindowCheckType	14
3.7.21	useTubeTexture	14
3.7.22	windowFailState	14
3.7.23	windowPassState	14
3.8	stateTableText	14
3.9	TrialTargets	14
3.10	TrialTubeParameters	15
3.11	visualLag	15
3.12	visualLagRandom	15

4	TrialData	16
4.1	analogData	16
4.2	Decoder	16
4.2.1	averagedSpikeBinHistory	16
4.2.2	decoderType	16
4.2.3	Parameters	16
4.2.4	rawDecode	16
4.2.5	rawSpikeBins	16
4.2.6	runDecoder	16
4.3	intertrialInterval	16
4.4	Marker	16
4.4.1	autoMonkey	16
4.4.2	errorCursor	17
4.4.3	frequency	17
4.4.4	rawPositions	17
4.4.5	rotation	17
4.4.6	rotationCenter	17
4.4.7	scaling	17
4.4.8	SyncParameters	18
4.4.9	trackingMarker	18
4.5	numberOfManualRewards	18
4.6	stateTransitions	18
4.7	TDT	18
4.7.1	bb_samplingRate	18
4.7.2	bb_waveforms	18
4.7.3	Block	18
4.7.4	brainArea	18
4.7.5	lfpCutoffFrequencies	18
4.7.6	lfpCutoffFrequenciesFs	18
4.7.7	lfpFs	19
4.7.8	lfpWaveforms	19
4.7.9	numberOfChannels	19
4.7.10	rawTDT	19
4.7.11	snippetInfo	19
4.7.12	snippetInfoCutoffFrequencies	19
4.7.13	snippetInfoCutoffFrequenciesFs	19
4.7.14	snippetWaveforms	19
4.7.15	sortName	19
4.7.16	tankPath	19
4.7.17	threshold	19
4.7.18	thresholdFs	19
4.7.19	thresholdStatic	19
4.7.20	thresholdStaticFs	20
4.7.21	timestampUnits	20
4.7.22	trialID	20

Data Structure (Version 2019.01.25)

1 Overview

Miscellaneous trial information.

1.1 date

Date of trial data collection (saved in numeric yyyyymmddHHMMSS format).

1.2 motionTracker

Which motion tracker was used for this trial (legacy from when we used wiimotes; if blank, assume phasespace was used).

1.3 parameterFile

File path to most recent .prm file loaded before this trial.

1.4 rigCodeVersion

Rig code executable version used to collect the data for this trial.

1.5 singleTrialTag

Trial was tagged (1) or not (0).

1.6 subjectName

Subject for this trial.

1.7 tag

Experiment block this trial was in.

1.8 translationDate

Date this trial was translated (saved in numeric yyyyymmddHHMMSS format).

1.9 trialName

State table name for this trial.

1.10 trialNumber

Trial number for this trial (saved as 'Trial[number]').

1.11 trialStatus

Status of this trial.

- 0 = Failed trial
- 1 = Successful trial
- 2 = Error occurred during trial

2 Definitions

Definitions of all possible parameter values for this trial.

2.1 `algTargTypeDef`

All data sources for analog targets.

2.2 `AllPossibleAnalogTargets`

Properties of all defined analog targets. Analog targets use voltage data (see [Section 4.1](#) for more information) to determine success.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one analog target.

2.2.1 `dataSource`

Enum value associated with data source (zero-indexed, listed in order in [Section 2.1](#)) for each analog target.

2.2.2 `dataSourceType`

Data source (defined in [Section 2.1](#)) for each analog target.

2.2.3 `linkedTargets`

Linked marker target for each analog target (only defined if `dataSourceType` is Gaze).

2.2.4 `names`

Name of each analog target.

2.2.5 `tolerances`

Number of time points analog window conditions can be not met before window conditions are determined to be not held for each analog target.

2.2.6 `windows`

Analog target tolerance windows.

For TTL and Touch analog targets, the first value is a voltage minimum threshold.

Gaze tolerance windows match the marker target tolerance window (see [Section 5](#)) of its linked target.

2.3 `allPossibleIntertrialIntervals`

List of all possible intertrial intervals (ITIs) (in "ms").

2.4 `AllPossibleIntervals`

List of all defined interval variables.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one interval.

2.4.1 modes

Value selected from list (0) or calculated from formula (1) (see [Section 2.21](#)).

2.4.2 names

Name of interval variable.

2.4.3 primary

Whether each interval can be used in state table steps (1) or is only valid for intermediate calculation (0).

2.4.4 values

Listed interval values. Only nonzero values are valid.

2.5 AllPossibleMarkers

Properties of all defined markers.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one marker.

2.5.1 channelNumber

Channel number for each marker (only valid for Phasespace markers).

2.5.2 dataSource

Enum value associated with data source for each marker

- 0 = Phasespace/Wiimotes
- 1 = Neural Decoded
- 2 = Auto-monkey
- 3 = Error
- 4 = Force

2.5.3 dataSourceType

Data source (Phasespace/Wiimotes, Neural Decoded, Auto-monkey, Error, Force) for each marker.

2.5.4 names

Name of each marker.

2.6 AllPossibleStateTables

Properties of all possible state tables.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one state table.

2.6.1 globalTimeout

Global timeout value for each state table (in "ms"). If global timeout is applied, then when the specified time has elapsed, the trial will jump to the state table step defined in globalTimeout-FailureStep.

2.6.2 globalTimeoutFailureStep

Step to jump to once global timeout period has elapsed for each state table.

2.6.3 names

Name of each state table.

2.6.4 probabilities

Probability for each state table (in percentage).

2.7 AllPossibleTargets

Properties of all possible marker targets.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one target.

2.7.1 formulas

Formula for calculating location of each marker target (valid only for Cartesian Formula mode marker targets).

2.7.2 modes

Mode of marker target location selection (see [Section 2.26](#)) for each marker target.

2.7.3 names

Name of each marker target.

2.7.4 primary

Whether each marker target can be referenced in state table steps (1) or is only valid for intermediate calculation (0).

2.7.5 selectionModes

Mode of selecting target locations from list (see [Section 2.27](#)).

2.7.6 shape

Shape of marker target window (see [Section 2.30](#)).

2.7.7 window

Tolerance window for each marker target (in mm).

For box-shaped (see [Section 2.30](#)) tolerance windows, these values represent x, y, and z tolerance values. For cylindrical tolerance windows, the numerical values represent radius and depth (the second value is NaN). For spherical tolerance windows, the numerical value represents radius (the second and third values are NaNs).

2.7.8 xValues

Lists of possible x position values (in mm).

2.7.9 yValues

Lists of possible y position values (in mm).

2.7.10 zValues

Lists of possible z position values (in mm).

2.8 AllPossibleTubeFilepaths

¹ File path to each tube.

2.9 AllPossibleTubeGroups

Group(s) assigned to each tube. Zero values are only valid if listed first.

2.10 AllPossibleTubeLengths

Default distance between endpoints of each tube.

2.11 AllPossibleTubeNames

Name of each tube.

2.12 AllPossibleTubePathLengths

Default path length for each tube

2.13 AllPossibleTubeRenderTypes

Graphical rendering type for each tube.

2.14 allPossibleVibrotactorDefinitions

Types of motors.

2.15 allPossibleVibrotactorDimensions

Possible directional modes of motor activation.

¹ Consider consolidating tube properties into struct.

2.16 allPossibleWaveformShapes

Types of tactor or pager motor waveform shapes.

2.17 analogChannelNames

Types of information saved in analogData (see [Section 4.1](#)) corresponding to each column.

2.18 autoMonkeyProfileDef

Types of auto-monkey trajectory profiles:

- Linear: Linear speed profile
- Bell-shaped: Bell-shaped speed profile
- Trajectory: User-defined times and positions

2.19 gazeCheckDef

Types of checks on eye position.

2.20 handCheckDef

Types of checks on marker targets (and analog targets of types other than TTL).

2.21 intervalModeDef

Modes of selecting or calculating interval values.

2.22 MarkerSpeedProfiles

Properties of all possible auto-monkey trajectory profiles.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one auto-monkey profile.

2.22.1 finalValues

Final cursor speed value (in mm/s if absolute speed) for each auto-monkey trajectory profile.

2.22.2 initialValues

Initial cursor speed value (in mm/s if absolute speed) for each auto-monkey trajectory profile.

2.22.3 names

Name of each auto-monkey trajectory profile

2.22.4 relativeSpeed

Whether speed is defined relative to start and end locations (1) or in absolute magnitude (0).

2.22.5 types

Type of trajectory for each auto-monkey trajectory profile (see [Section 2.18](#)).

2.23 motorModeDef

Activation modes for tactors or pager motors.

2.24 rawDecodeColumnNames

Types of information saved in rawDecode (see [Section 4.2.4](#)) corresponding to each column.

2.25 snippetInfoRowNames

Types of extracted TDT snippet information corresponding to each row.

2.26 targetModeDef

Modes of selecting or calculating target locations.

2.27 targetSelectionModeDef

Modes of selecting target locations from possible values (valid only for Cartesian List mode targets).

2.28 targetShapeDef

Graphical display shape options for marker targets.

2.29 targetTypeDef

Types of graphical display objects.

2.30 targetWindowShapeDef

Types of marker target tolerance window shapes.

2.31 ttlCheckDef

Types of checks on TTL voltage.

2.32 tubeCheckDef

Types of checks on tube window.

2.33 vrPhasespaceTransformation

Parameters for calculating VR space position from Phasespace coordinates.

2.33.1 rotation

Rotation about x, y, and z axes (in degrees).

2.33.2 scaling

Scaling of x, y, and z axes.

2.33.3 translation

Translation along x, y, and z axes (in mm).

3 Parameters

3.1 eyeTrackerCalibrationMatrix

Matrix mapping eye tracker analog voltage outputs to Phasespace coordinates. The first column represents x and y translation (in mm). The next two columns are gain multipliers for the left eye x and y voltage outputs (in mm/V). If binocular tracking was used, the last two columns are gain multipliers for the right eye x and y voltage outputs (in mm/V).

3.2 ForceParameters

CST error cursor parameters. ²

3.2.1 initialLambda

Initial lambda value for the dynamical system.

3.2.2 inputMethod

Data source used as input to dynamical system

- 0 = Force (i.e., pinch force)
- 1 = Decoded Cursor
- 2 = Phasespace

3.2.3 kicms

Intracortical microstimulation gain (no ICMS if 0).

3.2.4 kv

Visual gain (no visual feedback if 0).

3.2.5 kvib

Vibrotactile gain (no tactor feedback if 0).

3.2.6 refValue

Reference (target) value of the dynamical system. This value is subtracted from the raw input before input is passed through the dynamical system.

3.2.7 slowfastRate

Rates defining how lambda changes based on the error. The fast rate is used until the error surpasses a given threshold, at which point the slow rate is used.

3.2.8 systemName

Description of the dynamical system.

²Neither error threshold nor visual direction are saved to data structure. Additionally, kvib does not appear to be used online. If kv is 0, error cursor input appears to always be 0. In old vibrotactile CST data, is kv 1 but error cursor not displayed? Or was this written out in the dynamical system response code?

3.2.9 systemOrder

Order of the dynamical system (1 or 2).

3.2.10 visCursorOffsets

Offset applied to cursor for visual feedback (in mm).

3.3 LockZParameters

Parameters for overwriting Phasespace z coordinate.

3.3.1 lockZ

Overwrite Phasespace z coordinate (1) or not (0).

3.3.2 lockZCoordinate

Value to use in place of current Phasespace z coordinate (in mm).

3.4 markerJitter

Phasespace jitter parameters.

3.4.1 A

Gain multiplier for previous jitter value.

3.4.2 controlJitter

Use jittered cursor position for online task control (1) or not (0).

3.4.3 S

Standard deviation of jitter (in mm).

3.4.4 visibleJitter

Display cursor visually with jitter (1) or not (0).

3.5 MotorParameters

Properties of tactor and pager motor activation.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one motor.

3.5.1 amplitudeRanges

Minimum and maximum amplitudes for each motor for this trial (in V).

3.5.2 carrierFrequencyRanges

Minimum and maximum carrier frequencies for each motor (in Hz).

3.5.3 distanceRanges

Minimum and maximum distances for each motor (in mm).

3.5.4 dutyCycleRanges

Minimum and maximum duty cycle values for each motor (in percentage).

3.5.5 forceToZeroIfOutsideRange

Coerce value to 0 if outside range (1) or not (0) for each motor.

3.5.6 ignoreIfOutsideRange

Ignore value if outside range (1) or not (0) for each motor.

3.5.7 modulationFrequencyRanges

Minimum and maximum modulation frequencies for each motor (in Hz).

3.5.8 vibrotactorDimensions

Dimension for each motor (see [Section 2.15](#)).

3.5.9 vibrotactorsNumber

Tactor number for each motor.

3.5.10 vibrotactorTypes

Type (see [Section 2.14](#)) of each motor.

3.5.11 waveformShapes

Waveform shape (see [Section 2.16](#)) for each motor.

3.6 Pseudorandomness

3.6.1 pseudorandom

Pseudorandomize target presentation for this trial type (1) or not (0).

3.6.2 pseudorandomReset

Pseudorandomness was reset this trial (1) or not (0). ³

3.7 StateTable

Each element of this array represents the trial progression parameters for a single state of this trial. These array elements contain the following fields:

³*Pseudorandomness reset value is always reset to 0 before being sent back to Host.*

3.7.1 Analog

Properties of analog targets referenced in this state.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one analog target.

3.7.1.A. dataSource

Data source for each analog target.

3.7.1.B. linkedTarget

Name of linked target for each analog target.

3.7.1.C. linkedTargetIndex

Index of linked target, if valid, in TrialTargets array (see [Section 3.9](#)).

3.7.1.D. targetName

Name of analog target.

3.7.1.E. targetNameIndex

Index of analog target in AllPossibleAnalogTargets array (see [Section 2.2](#)).

3.7.1.F. targetWindowCheckType

Enum value of analog target window check type (see [Sections 2.20 and 2.31](#)).

3.7.1.G. window

Tolerance window for each analog target (see [Section 2.2.6](#)).

3.7.2 applyGlobalTimeout

Apply global timeout (1) or not (0) (see [Section 2.6.1](#)).

3.7.3 autoMonkeySpeedProfileName

Auto-monkey trajectory profile for this state (N/A if no auto-monkey trajectory profile selected).

3.7.4 Display

Graphical display properties for this state.

3.7.4.A. cursor

Display cursor (1) or not (0) for this state.

3.7.4.B. endTarget

Display endTarget (1) or not (0) for this state.

3.7.4.C. markers

Marker numbers to display and whether to freeze marker (1) or not (0) for this state.

3.7.4.D. startTarget

Display startTarget (1) or not (0) for this state.

3.7.4.E. tube

Display tube (1) or not (0) for this state.

3.7.5 Gaze

Gaze fixation target for this state.

3.7.5.A. targetName

Name of linked gaze target.

3.7.5.B. targetNameIndex

Index of linked target, if valid, in TrialTargets array.

3.7.5.C. targetWindowCheckType

Enum value of gaze target window check type (see [Section 2.20](#)).

3.7.5.D. window

Tolerance window for the linked gaze target (see [Section 5](#)).

3.7.6 Hand

Marker target properties for this state.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one marker target.

3.7.6.A. targetMarker

Name of marker used for window checks for this target.

3.7.6.B. targetName

Name of marker target.

3.7.6.C. targetNameIndex

Index of marker target in TrialTargets array.

3.7.6.D. targetWindowCheckType

Enum value of marker target window check type (see [Section 2.20](#)).

3.7.6.E. window

Tolerance window for each marker target.

Column format:

- x coordinate of target center (in mm)
- y coordinate of target center (in mm)
- z coordinate of target center (in mm)
- x window length (see [Section 2.7.7](#)) (in mm)
- y window length (in mm)
- z window length (in mm)

3.7.7 Interval

Interval used for this state.

3.7.7.A. length

Value of interval (in "ms").

3.7.7.B. name

Name of interval.

3.7.8 Motor

Motor control parameters for this state.

3.7.8.A. mode

Mode of motor activation (see [Section 2.23](#)).

3.7.8.B. target

Location of motor target (x, y, z coordinates, in mm).

Note: If no motor target defined, location will be reported as [-1 -1 -1].

3.7.8.C. **targetName**

Name of motor target.

3.7.8.D. **targetNameIndex**

Index of motor target in TrialTargets array.

3.7.9 **passTable**

Define the trial as successful if state is reached (1) or not (0).

3.7.10 **resetGlobalTimeout**

Reset global timeout counter (1) or not (0) (see [Section 2.6.1](#)) for this state.

3.7.11 **stateName**

Name of this state.

3.7.12 **StateTargets**

Displayed target properties for this state.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one displayed target.

3.7.12.A. **color**

Color of displayed target.

3.7.12.B. **location**

Location of displayed target (x, y, z coordinates, in mm).

3.7.12.C. **nameIndex**

Index of displayed target in TrialTargets array.

3.7.12.D. **names**

Name of displayed target.

3.7.12.E. **shape**

Shape of displayed target (see [Section 2.28](#)).

3.7.12.F. **size**

Size of displayed target (in mm).

3.7.12.G. **type**

Graphical target type (see [Section 2.29](#)) of displayed target.

3.7.13 **togglePhotodiodeTarget**

Toggle photodiode target (1) or not (0) for this state.

3.7.14 **TTL**

TTL output check type and output pattern for this state. ⁴

⁴ *Output pattern names and values not saved to data structure.*

3.7.15 tubeColor

Color of tube for this state (only valid if tube is displayed).

3.7.16 tubeMarker

Marker for tube control for this state (only valid if tube window is checked).

3.7.17 tubeSection

Section of tube for control and display for this state (if 0, the whole tube is used for control and display).

3.7.18 tubeTexture

Texturing applied to tube for this state.

3.7.19 tubeWindow

Tube tolerance radius and depth for this state (in mm).

3.7.20 tubeWindowCheckType

Tube window check type for this state (see [Section 2.32](#)).

3.7.21 useTubeTexture

Use tube texture (1) or not (0) for this state.

3.7.22 windowFailState

State to jump to if this state fails (0 represents trial end).

3.7.23 windowPassState

State(s) to jump to if this state succeeds (0 represents trial end). ⁵

3.8 stateTableText

All state table cell values for this trial.

3.9 TrialTargets

Properties of all displayed targets for this trial.

Note: Each of the following fields contains an array, where each vector element (or matrix row) corresponds to one displayed target.

3.9.0.A. color

Color of each target.

3.9.0.B. names

Name of each target.

3.9.0.C. size

Size of each target.

⁵State table target logic (AND vs OR) not saved in data structure.

3.9.0.D. statesUsed

States where each target was used.

3.9.0.E. window

Tolerance window of each displayed target (see [Section 5](#)).

3.10 TrialTubeParameters

3.10.0.A. displayRadii

Graphical radius of displayed tube (in mm).

3.10.0.B. numberSections

Total number of sections of tube.

3.10.0.C. pseudorandomizeTubes

Pseudorandomize tube presentation (1) or not (0).

3.10.0.D. pseudorandomTubeTarget

Name of target linked to tube pseudorandomness.

3.10.0.E. rotation

Rotation of target trajectory (in degrees).

3.10.0.F. size

Size of tube (valid only for very old data).

3.10.0.G. trajectory

Downsampled tube trajectory coordinates (in mm).

3.10.0.H. tubeName

Name of tube.

3.10.0.I. tubeType

Type of tube (valid only for very old data).

3.10.0.J. window

Radius and depth of tube control window (in mm).

3.11 visualLag

Visual lag applied to phasespace data (in "ms").

3.12 visualLagRandom

Randomize visual lag (1) or not (0).

4 TrialData

Data used for task control for this trial.

4.1 analogData

Voltage data from analog input channels. Each row represents one frame.

The "actual time" column collected in newer executables contains the number of 500 us periods elapsed from the end of the last frame. Thus the actual time in ms can be calculated by multiplying this column by 0.5.

4.2 Decoder

4.2.1 averagedSpikeBinHistory

Average spike count history for this trial.

4.2.2 decoderType

Type of decoder (e.g. Kalman (for data collected with older executables), LDS, GPFA) for this trial.

4.2.3 Parameters

Parameters of decoder for this trial.

4.2.4 rawDecode

Raw decode information for this trial. Each row represents one decoder prediction. Column information can be found in [Section 2.24](#).

4.2.5 rawSpikeBins

Raw spike bin information for this trial. Raw spike bins used online for decoder prediction.

4.2.6 runDecoder

Run decoder (1) or not (0) for this trial.

4.3 intertrialInterval

Intertrial interval for this trial (in "ms").

4.4 Marker

4.4.1 autoMonkey

Auto-monkey cursor data for this trial. The auto-monkey follows a user-defined trajectory.

4.4.1.A. position

Auto-monkey cursor positions (in mm).

4.4.1.B. time

Auto-monkey cursor time (in "ms").

4.4.2 errorCursor

CST error cursor data for this trial. The error cursor is defined by a dynamical system.

Column format:

- Error cursor x coordinate (in mm)
- Error cursor y coordinate (in mm)
- Error cursor z coordinate (in mm)
- Time (in "ms")
- Lambda value (in rad/s)
- Dynamical system input value (in mm)

4.4.3 frequency

Phasespace frequency for this trial.

4.4.4 rawPositions

Phasespace position data for this trial.

Column format:

- Phasespace marker ID
- Phasespace raw x coordinate (in mm)
- Phasespace raw y coordinate (in mm)
- Phasespace raw z coordinate (in mm)
- Phasespace frame number
- Phasespace synchronized raw time (in ms, calculated from frame number)
- Phasespace synchronized raw time + visual lag (in ms)
- Phasespace x jitter (in mm)
- Phasespace y jitter (in mm)
- Phasespace z jitter (in mm)
- Phasespace jittered x coordinate (in mm)
- Phasespace jittered y coordinate (in mm)
- Phasespace jittered z coordinate (in mm)

4.4.5 rotation

Phasespace rotation value (in degrees).

4.4.6 rotationCenter

Coordinates of phasespace rotation origin (in mm).

4.4.7 scaling

Scaling of phasespace data.

4.4.8 SyncParameters

4.4.8.A. phasespaceFrame

Phasespace frame number value at start of sync pulse.

4.4.8.B. pulseLength

Length of sync pulse from start time (in "ms").

4.4.8.C. pulseLogic

Digital line output value of sync pulse (1 = 5V, 0 = 0V). Outside of sync pulse, digital output line is set to opposite binary value.

4.4.8.D. startTime

Start time of sync pulse (in "ms").

4.4.9 trackingMarker

Phasespace marker number used for movement tracking.

4.5 numberOfManualRewards

Number of manual rewards given via pushbutton.

4.6 stateTransitions

State transition data. The first row is state number. The second row is the frame number at which this state began. The last state (-1) represents the end of the trial.

4.7 TDT

Neural data saved for this trial.

4.7.1 bb_samplingRate

Sampling rate for broadband waveform for this trial (in Hz).

4.7.2 bb_waveforms

Broadband waveform data for this trial (in V).

4.7.3 Block

Name of original block in TDT tank where neural data was saved for this trial.

4.7.4 brainArea

Brain area where this trial's subject's array was implanted.

4.7.5 lfpCutoffFrequencies

Local field potential (LFP) filter frequencies for this trial (in Hz).

4.7.6 lfpCutoffFrequenciesFs

Sampling rate for LFP filter frequencies for this trial (in Hz).

4.7.7 lfpFs

Sampling rate for LFP data for this trial (in Hz).

4.7.8 lfpWaveforms

LFP waveform data for this trial (in V).

4.7.9 numberOfChannels

Number of channels input to the PZ2 preamplifier for this trial.

4.7.10 rawTDT

Vector of UDP packets of neural data from the RZ2 for this trial.

Each packet contains a header, RT time, TDT sample number, and 32-bit words representing spike counts of multiple channels.

4.7.11 snippetInfo

Selected snippet event information for this trial (see [Section 2.25](#)).

4.7.12 snippetInfoCutoffFrequencies

Neural filter cutoff frequencies for this trial (in Hz).

4.7.13 snippetInfoCutoffFrequenciesFs

Sampling rate for neural filter cutoff frequencies for this trial (in Hz).

4.7.14 snippetWaveforms

Snippet waveform data for this trial (in V).

4.7.15 sortName

Name of sort data for this trial.

4.7.16 tankPath

Path to tank directory for neural data stored for this trial.

4.7.17 threshold

Threshold values for this trial (in V).

4.7.18 thresholdFs

Sampling rate for threshold values for this trial (in V).

4.7.19 thresholdStatic

Manual (1) or automatic (0) threshold values for this trial.

4.7.20 thresholdStaticFs

Sampling rate of manual threshold status for this trial (in Hz).

4.7.21 timestampUnits

Neural event timestamp units for this trial. Millisecond (ms) or microsecond (us) precision can be used.

4.7.22 trialID

Trial ID value sent to TDT to mark beginning of this trial.